
EMORY UNIVERSITY
Environmental Health and Safety Office
Research Administration

1762 Clifton Road NE, Suite 1200
Atlanta, Georgia 30322
(404) 727-5922
FAX: (404) 727-9778


PROCEDURE FOR ACCESS TO INDIVIDUAL DOSE REPORT



- Individual Dose Report allows the client to setup a generic Username/Password that can be given to all participants and allow them to use myLDR to look up their dose history.
- The participant will use the serial number located on the back of their dosimeter to retrieve information. Only the serial number located on the back of the dosimeter will work because it includes an extra character. This is for added security.
- When the participant accesses the page there is no participant information listed for added security as well.



- This summary screen allows the participant to view their most recent quarter's results, the last 2 year's totals, and lifetime.
- If participant would like to see more detailed information they can click the View Details button.



History Results

The doses are displayed in mrem

Dose Results

Total Records: 4

Doses as of 2013/03/13 20:12 CST

Dose Period	Total DDE	Total LDE	Total SDE	Beta	Total Neutron	Extremity
Q32012	M	M	M			
2012	11	11	10			
2011	234	234	225			
Lifetime	888	888	859			

Page 1 of 1 Prev Next

Export Result to : Excel Export View Details

History Details


Dosimeter Dose

Total Records: 24

Doses as of 2013/03/13 20:12 CST

Account	Subaccount Name	Subaccount Code	Begin Wear Date	End Wear Date	Dosimeter Type	Dosimeter Location	Total DDE	Total LDE	Total SDE	Beta	Total Neutron	Fast Neutron	Thermal Neutron	Extremi
182051	RAD	RAD	2013/01/01	2013/01/31	Pa	Chest								
182051	RAD	RAD	2012/12/01	2012/12/31	Pa	Chest								
182051	RAD	RAD	2012/11/01	2012/11/30	Pa	Chest								
182051	RAD	RAD	2012/10/01	2012/10/31	Pa	Chest								
182051	RAD	RAD	2012/09/01	2012/09/30	Pa	Chest	M	M	M					
182051	RAD	RAD	2012/08/01	2012/08/31	Pa	Chest	M	M	M					
182051	RAD	RAD	2012/06/01	2012/06/30	Pa	Chest	M	M	M					
182051	RAD	RAD	2012/05/01	2012/05/31	Pa	Chest	M	M	M					

- The View Details button allows the participant to see their exposure month by month for the last 2 years.
- Fetal doses are not included in this history because they do not carry assigned doses. Only the mother's doses are present.



A	B	C	D	E	F	G	H	I	J
Account	Subaccount Name	Subaccount Code	Begin Wear Date	End Wear Date	Dosimeter Type	Dosimeter Location	Total DDE	Total LDE	Total SDE
182051	RAD	RAD	2013/01/01	2013/01/31	Pa	Chest			
182051	RAD	RAD	2012/12/01	2012/12/31	Pa	Chest			
182051	RAD	RAD	2012/11/01	2012/11/30	Pa	Chest			
182051	RAD	RAD	2012/10/01	2012/10/31	Pa	Chest			
182051	RAD	RAD	2012/09/01	2012/09/30	Pa	Chest	M	M	M
182051	RAD	RAD	2012/08/01	2012/08/31	Pa	Chest	M	M	M
182051	RAD	RAD	2012/06/01	2012/06/30	Pa	Chest	M	M	M
182051	RAD	RAD	2012/05/01	2012/05/31	Pa	Chest	M	M	M
182051	RAD	RAD	2012/04/01	2012/04/30	Pa	Chest	4	4	4
182051	RAD	RAD	2012/03/01	2012/03/31	Pa	Chest	4	4	4
182051	RAD	RAD	2012/02/01	2012/02/29	Pa	Chest	3	3	2
182051	RAD	RAD	2012/01/01	2012/01/31	Pa	Chest	M	M	M
182051	RAD	RAD	2011/12/01	2011/12/31	Pa	Chest	M	M	M
182051	RAD	RAD	2011/11/01	2011/11/30	Pa	Chest	29	29	29
182051	RAD	RAD	2011/10/01	2011/10/31	Pa	Chest	M	M	M
182051	RAD	RAD	2011/09/01	2011/09/30	Pa	Chest	29	29	27
182051	RAD	RAD	2011/08/01	2011/08/31	Pa	Chest	31	31	31
182051	RAD	RAD	2011/07/01	2011/07/31	Pa	Chest	M	M	M
182051	RAD	RAD	2011/06/01	2011/06/30	Pa	Chest	M	M	M
182051	RAD	RAD	2011/05/01	2011/05/31	Pa	Chest	M	M	M
182051	RAD	RAD	2011/04/01	2011/04/30	Pa	Chest	33	33	32
182051	RAD	RAD	2011/03/01	2011/03/31	Pa	Chest	35	35	33
182051	RAD	RAD	2011/02/01	2011/02/28	Pa	Chest	35	35	33
182051	RAD	RAD	2011/01/01	2011/01/31	Pa	Chest	42	42	40

- All dose information can be exported.